



AVISCERA BIOSCIENCE

Rabbit Anti-Human ATPase Inhibitory factor 1 (ATPIF1) Purified IgG

Product Information

Code	A00598-01-100
Name	Anti ATPIF1 (Human)
Clone No.	polyclonal
Lot No.	
Size	100 µg
Species	Human
Host	Rabbit
Immunogen	human ATPIF1 rec
Ab Type	IgG
Purification	Protein A
Formulation	lyophilized
Carry	Carry free
Storage	-20 ° C
Specificity	human
Reconstitution	PBS, 100 µl
Application	IHC WB

Preparation

This antibody was produced from a rabbit immunized with highly purified recombinant human ATPIF1. The antibody was purified by Protein A.

Formulation

100 µg of rabbit anti-human ATPIF1 IgG in PBS, pH 7.4. Carry Free.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize human mitochondrial ATPase Inhibitor 1 on ELISA, Western Blot and Immunohistochemistry.

Applications

Western blot - This antibody can be used at 2~4 µg/mL with the appropriate secondary reagents to detect recombinant human ATPIF1 under reduce condition.

Immunohistochemistry -That Antibody can be used at 4 µg/mL with the appropriate secondary antibody to detect ATPIF1 on paraffin-embedded human heart tissue and colon cancer tissue (ABC).

ELISA-This antibody can be used as Capture Antibody at 1-2 µg/ml combine with Detection Antibody A00598-02-50B to measure recombinant human ATPIF1 on ELISA.

Optimal dilutions should be determined by each laboratory for each application.

ORDER INFORMATION

AVISCERA BIOSCIENCE INC.

2348 Walsh Ave. Suite C

Santa Clara, CA 95051

Tel: (408) 982 0300

Fax: (408) 982 0301

Email:

Sales@AvisceraBioscience.com

www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.