



Rabbit Anti-Mouse Cerebellin-1 IgG biotinylated

Product Information

Code	A00361-03-50B
Name	CBLN1 (M) Pab
Clone No.	N/A
Lot No.	
Size	50 µg
Species	Mouse
Host	Rabbit
Immunogen	Cerebellin-1 (Mouse), mature form
Ab Type	IgG
Purification	Protein A
Formulation	Lyophilized Form without Preservatives
Carry	Free
Storage	-20 ° C
Specificity	M, R, H
Reconstitution	100 µl
Application	IHC ELISA

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Sales@AvisceraBioscience.com
www.AvisceraBioscience.com

Preparation

This antibody was produced from a rabbit immunized with purified, *E. coli*-derived, recombinant mature form of mouse Cerebellin-1. That IgG was purified by Protein A affinity and conjugated with water soluble biotin.

Formulation

50µg of purified anti mouse cerebellin-1 IgG in 50 µl of TBS contain 0.1% BSA without preservatives was lyophilized.

Reconstitution

Add 50 µl of PBS to the vial to prepare antibody stock solution at 0100 µg/mL. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for up six months without detectable loss of activity.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize mature form of mouse Cerebellin-1 on indirect ELISAs.

Applications

Indirect ELISA - This antibody can be used at 1: 8000 (0.125 µg/mL) with the appropriate secondary reagents to detect Cerebellin-1 (Human) on indirect ELISA and immunohistochemistry.

Immunohistochemistry-That Antibody can be used at 2-4 µg /mL with the appropriate secondary antibody to detect Cerebellin-1 in the paraffin embedded mouse, rat and human brain tissues (ABC).

Optimal dilutions should be determined by each laboratory for each application.