



Anti Human Follistatin-Related Protein 1 (FSTL1) Monoclonal Antibody

Product Information

Code	A00347-24-100
Name	Anti FSTL1 (Human) Mab
Clone No.	3D10A1A11
Lot No.	
Size	100 µg
Species	Human
Host	Mouse
Immunogen	Human FSTL1, rec.
Ab Type	IgG
Purification	Protein G
Formulation	Lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human
Reconstitution	100 µl
Application	ELISA

ORDER INFORMATION

AVISCERA BIOSCIENCE, Inc

2348 Walsh Ave, Suite C,
Santa Clara, CA 95051
USA

Tel: +001 408 982 0300

Fax: +001 408 982 0301

Info@Aviscerabioscience.com

www.Aviscerabioscience.com

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified recombinant human FSTL1. That antibody was purified by Protein G affinity.

Formulation

100 µg of Anti Human FSTL1 monoclonal antibody in 100 µl of PBS lyophilized form.

Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (1000µg/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for 12 months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize recombinant human FSTL1 on indirect ELISA.

Applications

Indirect ELISA - This antibody can be used at 0.031µg/ml to detect human FSTL1 on indirectly ELISA.

ELISA - This antibody can be used as capture antibody at 1-2 µg/ml paired with detection antibody rabbit Anti FSTL1 (Human) IgG (A00347-02-100) to detect human FSTL1 (00347-02-100) on ELISA.

Optimal dilutions should be determined by each laboratory for each application.

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.