



AVISCERA BIOSCIENCE

Anti-Rat Fetuin A IgG

Product Information

Code	A00173-02-100
Name	Anti Rat Fetuin A Antibody
Clone No.	N/A
Lot No.	
Size	100 µg
Species	Rat
Host	Rabbit
Immunogen	Rat Fetuin A rec.
Ab Type	IgG
Purification	Protein A affinity
Formulation	lyophilized Form without preservatives free
Carry	
Storage	-20~-70 ° C
Specificity	Rat
Reconstitution	PBS, 100 µl
Application	IHC, ELISA

PREPARATION

This antibody was produced from a rabbit immunized with purified, *E. coli*-derived, recombinant rat Fetuin A. That IgG was purified by Protein A affinity

FORMULATION

100 µg of Rabbit Anti Rat Fetuin A purified IgG in 100 µl of PBS without preservatives was lyophilized.

RECONSTITUTION

Add 100 µl of PBS to the vial to prepare antibody stock solution at 100 µg/100 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for up six months without detectable loss of activity.

STORAGE

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ~ -70 ° C for 12 months. **Avoid repeated freeze-thaw cycles.**

SPECIFICITY

This antibody has been selected for its ability to recognize rat Fetuin A in indirect ELISA and immunohistochemistry.

APPLICATIONS

Indirect ELISA - This antibody can be used at 1:8000 (0.125 µg/ml) to detect recombinant rat Fetuin A on indirect ELISA.

Immunohistochemistry- That antibody can be used at 2-4µg/ml to detect Fetuin A in rat liver tissue slides (ABC).

Optimal dilutions should be determined by each laboratory for each application.

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

ORDER INFORMATION

AVISCERA BIOSCIENCE, INC
2348 Walsh Ave., Suite C
Santa Clara, CA 95051
Tel: (408) 982 0300
Info@Aviscerabioscience.com