



AVISCERA BIOSCIENCE

Rabbit Anti-Human Pro-GRP Antibody

Product Information

Code	A00320-01-100
Name	Anti Pro-GRP (Human) IgG
Clone No.	Polyclonal
Lot No.	
Size	100 µg
Species	Human
Host	Rabbit
Immunogen	Pro-GRP Rec
Ab Type	IgG
Purification	Affinity
Formulation	Dry
Carry	NO
Storage	-20 ° C
Specificity	Human
Reconstitution	PBS, 100 µl
Application	IHC
	E

Preparation

This antibody was produced from a rabbit immunized with Human Pro-GRP (54-148) and purified with Protein A affinity.

Formulation

Lyophilized 100 µg of rabbit anti-human Pro-GRP in 100 µl of PBS without any carry proteins.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize human Pro-GRP in indirect ELISA as well as immunohistochemistry.

Applications

Indirect ELISA - This antibody can be used at 0.125 µg/ml with the appropriate secondary reagents to detect human Pro-GRP (54-148) recombinant on indirect ELISA .

Immunohistochemistry -That Antibody can be used at 2-4 µg/mL with the appropriate secondary antibody to detect paraffin-embedded human lung cancer (ABC).

Optimal dilutions should be determined by each laboratory for each application.

ORDER INFORMATION

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Sales@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.