




### Detecting His Tagged Proteins in 45 minutes

### His Tag ELISA Kit SK00779-01



#### Feature:

- **Save Time & Cost:** Sample incubation for 30 min plus color developing for 10-15 min. \$ 3.75 per sample (80 samples individual test);
- Monoclonal Antibody Validated by WB
- One Step Assay in 45 min
- Calibrated by His Tagged Protein
- High Sensitivity Screening His Tagged Proteins
- 12 x 8-well Strips Easy for Partial Use
- Fast Screening His Tagged Expression and Purification
- Specificity to C-Terminus or N-Terminus His Tagged

#### His Tag ELISA Kit

Code No.: SK00779-01

Size: 96 T | 192 T

Price: \$298.00 USD | \$560.00 USD

Standard range: 1.6 ~ 1000 ng/ml

Sensitivity: 200 pg/ml

Specificity: 6 x, 5 x or 4 x His Tag on N-Terminal or C-Terminal


Calibration: His Tagged Protein Rec (11.7KD)

Sample Type: recombinant His Tagged proteins, superant of animal free cell cultures

Intra-CV: 4-8%

Inter-CV: 8-12%

#### His Tag ELISA Kit SK00779-01




4-PL Fit:  $y = (A - D) / (1 + (x/C)^B) + D$

|                                  | A    | B    | C    | D      | R <sup>2</sup> |
|----------------------------------|------|------|------|--------|----------------|
| STD#1 (Standards: Conc vs AvgOD) | 1.86 | 1.32 | 15.7 | 0.0598 | 0.999          |

Weighting: Fixed

## His Tag ELISA Kit - two Incubations, one hour assay



### His Tag ELISA Kit

Code No.: SK00779-06

Size: 96 T      480 T

Price: \$298.00 USD      \$1245.00 USD

Standard range: 1.6 ~ 200 ng/ml

Sensitivity: 200 pg/ml

Specificity: 6 x, 5 x or 4 x His Tag on N-Terminal or C-Terminal

Calibration: His Tagged Protein Rec (11.7KD)


Sample Type: recombinant His Tagged proteins, supernatant of animal free cell cultures

Intra-CV: 4-8%

Inter-CV: 8-12%

Data Sheet: PDF

**Pre-incubation of Standard or Samples with Anti His Tag Antibody in each well of PCR Micro-Plate 779-06-00**



## Custom Service: His Tagged Protein ImmunoAssay

### His Tagged Protein Assay

Code No.: CS00779-09

Size: 40 Samples      80 Samples

Price: \$1200.00 USD      \$1900.00 USD

Sample Type: purified recombinant His Tagged proteins, supernatant of animal free cell cultures

## His Tag ELISA Kit & Antibody

| Product Name                          | Applications                         | Code          | Size    | Price (\$) |
|---------------------------------------|--------------------------------------|---------------|---------|------------|
| His Tag ELISA Kit (96T)               | His tagged proteins or cell cultures | SK00779-01    | 96 T    | 298.00     |
| His Tag ELISA Kit (192T)              | His tagged proteins or cell cultures | SK00779-01F   | 192 T   | 560.00     |
| His Tag ELISA Kit (48T)               | His tagged proteins or cell cultures | SK00779-01H   | 48 T    | 249.00     |
| His Tag ELISA Kit (96T)               | His tagged proteins or cell cultures | SK00779-06    | 96 T    | 298.00     |
| His Tag ELISA Kit (480T)              | His tagged proteins or cell cultures | SK00779-06F   | 480 T   | 1245.00    |
| His Tag Western Blot Kit              | His tagged proteins or cell cultures | WBK00779-01   | each    | 650.00     |
| His Tag Antibody for Western Blot     | His tagged proteins or cell cultures | A00779-02-100 | 100 µL  | 369.00     |
| His Tag Antibody HRP for Western Blot | His tagged proteins or cell cultures | A00779-03-100 | 100 µL  | 369.00     |
| His Tagged Protein (11.7 KD) Rec      | Positive control                     | 779-01-04     | 1µg     | 99.00      |
| His Tagged Protein (11.7 KD) Rec      | Positive control                     | 779-01-04F    | 5 x 1µg | 460.00     |