



AVISCERA BIOSCIENCE

Anti Human IL-1F9 Monoclonal Antibody

Product Information

Code	A00739-06-100
Name	Anti IL-1F9 Mab
Clone No.	118-5-2G9
Lot No.	
Size	100 µg
Species	Human
Host	Mouse
Immunogen	Human IL-1F9, rec.
Ab Type	IgG
Purification	Protein G Lyophilized
Formulation	Form without preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human
Reconstitution	100 µl
Application	ELISA

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, E. coli-derived, recombinant human IL-1F9, His Tag on N-Terminus. That antibody was purified by Protein G affinity.

Formulation

100 µg of Anti Human IL-1F9 monoclonal antibody in 100 µl of PBS lyophilized form.

Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (100µg/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for 12 months without detectable loss of activity.

Avoid repeated freeze-thaw cycles.

Specificity

This antibody has been selected for its ability to recognize recombinant human IL-1F9 in indirect ELISA.

Applications

Indirect ELISA - This antibody can be used at 1:4000 (0.25 µg/ml) to detect human IL-1F9 on indirectly ELISA.

Optimal dilutions should be determined by each laboratory for each application.

ORDER INFORMATION

AVISCERA BIOSCIENCE, Inc

2348 Walsh Ave, Suite C,
Santa Clara, CA 95051

Tel: +001 408 982 0300

Fax: +001 408 982 0301

Info@Aviscerabioscience.com

www.Aviscerabioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.