



AVISCERA BIOSCIENCE

Human Soluble Interleukin 22 Binding Protein Recombinant

Product Information

Code	00527-01-50
Name	Human IL-22BP rec.
Lot No.	
Size	50 µg
Species	Human
Sequence	mature form
Protein ID	NP_851826
Gene ID	NM_181309
MW	26 KD
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS-PAGE gel PBS lyophilized
Formulation	Form without preservatives
Carry	free
Storage	-20 ° C ~ -70° C
Reconstitution	100 µl
Application	ELISA Cell biology

Description

A DNA sequence encoding the mature form of human soluble Interleukin-22 Binding Protein Isoform 2 (Thr²²-Pro²³¹) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

Formulation

Lyophilized 50 µg human soluble IL-22BP 50 µl of TBS, Carry free.

Reconstitution & Storage

Add 100 µl TBS to the vial to prepare a working stock solution at 50 µg/100 µl. Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: human soluble Interleukin-22 Binding Protein Isoform 2
(Thr²²-Pro²³¹)

AVISCERA BIOSCIENCE
2348 Walsh Ave., Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Fax: (408) 982 0301
Info@Aviscerabioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.