



Rabbit Anti- Human IL-33 IgG

Product Information

Code	A00277-01-100
Name	Human IL33 Pab
Clone No.	N/A
Lot No.	
Size	100 µg
Species	Human
Host	Rabbit
Immunogen	Human IL33 Rec.
Ab Type	IgG
Purification	Protein A
Formulation	Lyophilized Form without Preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human
Reconstitution	PBS, 100 µl
Application	ELISA

Preparation

This antibody was produced from a rabbit immunized with recombinant human IL-33. That IgG was purified by Protein A affinity.

Formulation

100 µg of purified IgG in 100 µl of PBS without preservatives was lyophilized.

Reconstitution

Add 100 µl of PBS to the vial to prepare antibody stock solution at 100 µg/100 µl. Store reconstituted antibody at 2 to 8 ° C for up to a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for up to six months without detectable loss of activity.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize mature Human IL-33 in direct ELISAs.

Applications

Indirect ELISA - This antibody can be used at 1: 8000 (0.125 µg /mL) with the appropriate secondary reagents to detect human IL33 on indirect ELISA.

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Info@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

Carry	free
Storage	-20 ° C
Specificity	Human, Rat