



Anti Human Interleukin 33 (IL-33) Monoclonal Antibody

Product Information

Code	A00277-14-100
Name	Anti IL-33 (H) MAB
Clone No.	1E8C8A11
Lot No.	
Size	100 µg
Species	Human
Host	Mouse
Immunogen	Human IL-33, rec.
Ab Type	IgG
Purification	Protein G PBS
Formulation	Lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human
Reconstitution	100 µl
Application	ELISA

Preparation

This antibody was produced from an hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, E. coli-derived, recombinant human IL-33. That antibody was purified by Protein G affinity.

Formulation

100 µg of Anti IL-33 (Human) Monoclonal Antibody in 100 µl of PBS lyophilized form.

Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (1000µg/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for 12 months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize recombinant human IL-33 on indirect ELISA and ELISA.

Applications

Indirect ELISA- Indirect ELISA - This antibody can be used at 1: 20000~40000 (0.025 ~ 0.05 µg /mL) with the appropriate secondary reagents to detect recombinant human IL-33 on indirect ELISA.

ELISA Assay - This antibody can be used as a capture antibody in a human IL-33 sandwich immunoassay in combination with the Rabbit Anti human IL-33 detection antibody (Code No.: A00277-01-100) and recombinant human IL-33 (Code No.: 00277-01-100) as the standard. The suggested concentration range for this capture antibody is 0.5-1 µg/mL in PBS.

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Sales@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.