



AVISCERA BIOSCIENCE

Anti-Human VEGF-B Monoclonal Antibody

Product Information

Code	A00122-06-100
Name	Human VEGF-B Mab
Clone No.	2E7B8
Lot No.	VEGFB006
Size	100 µg
Species	Human
Host	Mice
Immunogen	VEGF-B (H), rec.
Ab Type	IgG
Purification	<i>sequential precipitation</i>
Formulation	lyophilized Form without preservatives free
Carry	free
Storage	-20 ° C
Specificity	Human only
Reconstitution	PBS, 100 µl
Application	ELISA

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, *E. coli*-derived, recombinant human VEGF-B.

That antibody was purified by sequential precipitation caprylic acid and ammonium sulphate.

Formulation

100µg of Anti Human VEGF-B IgG in 100 µl of PBS without preservatives was lyophilized.

Reconstitution

Add 100 µl of PBS to the vial to prepare antibody stock solution at 100 µl.

Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for up six months without detectable loss of activity.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize human VEGF-B on direct ELISAs.

Applications

Indirect ELISA - This antibody can be used at 1: 128000 with the appropriate secondary reagents to detect human VEGF-B on indirect ELISA.

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE
2348 Walsh Ave., Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Info@AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

Biomarker Technology Solutions