



# AVISCERA BIOSCIENCE

## Rabbit Anti- von Willebrand Factor (vWF) (Human) IgG

### Product Information

Code	A00552-01-100
Name	Rabbit Anti vWF (H) IgG
Clone No.	Polyclonal
Lot No.	20110322
Size	100 µg
Species	Human
Host	Rabbit
Immunogen	Human vWF
Ab Type	IgG
Purification	Protein A
Formulation	Lyophilized Form without Preservatives
Carry	Free
Storage	-20 ° C
Specificity	human
Reconstitution	PBS, 100 µl
Application	IHC  ELISA

### PREPARATION

This antibody was produced from a rabbit immunized with purified, *E. coli*-derived, recombinant mature form of vWF (Human). That IgG was purified by Protein A affinity

### FORMULATION

100 µg of Rabbit Anti vWF (Human) purified IgG in 100 µl of PBS without preservatives was lyophilized.

### RECONSTITUTION

Add 100 µl of PBS to the vial to prepare antibody stock solution at 100 µg/100 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted ( by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for up six months without detectable loss of activity.

### STORAGE

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ~ -70 ° C for 12 months. **Avoid repeated freeze-thaw cycles.**

### SPECIFICITY

This antibody has been selected for its ability to recognize human vWF in indirect ELISAS , Western Blot and immunohistochemistry.

### APPLICATIONS

**Indirect ELISA** - This antibody can be used at 1:4000 (0.25 µg/ml) to detect recombinant human vWF on indirect ELISA.

**Western blot** - This antibody can be used at 5~10 µg/mL to detect recombinant human vWF under reduce condition.

**Immunohistochemistry-** That antibody can be used at 2~4 µg/ml to detect vWF in human lung and heart issue slide.

*Optimal dilutions should be determined by each laboratory for each application.*

AVISCERA BIOSCIENCE  
2348 Walsh Ave., Suite C  
Santa Clara, CA 95051  
USA  
Tel: (408) 982 0300  
Info@AvisceraBioscience.com

**THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.**