



AVISCERA BIOSCIENCE

Anti Human Secreted Frizzled-Related Protein 3 (sFRP3) IgG

Product Information

Code	A00308-06-100
Name	Human sFRP3 Pab
Clone No.	N/A
Lot No.	20111356
Size	100 µg
Species	Human
Host	Rabbit
Immunogen	sFRP3 (Human) Rec.
Ab Type	IgG
Purification	Protein A
Formulation	500µl
Carry	0.1%gelatin
Storage	-20 ° C
Specificity	H, M, R
Reconstitution	500 µl
Application	IHC, WB, IFC

Preparation

This antibody was produced from a rabbit immunized with recombinant human sFRP3 (156-325). That IgG was purified by Protein A affinity.

Formulation

100 µg of purified Rabbit Anti sFRP-3 IgG in 500 µl of PBS contains 0.1% gelatin.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize human, mouse and rat sFRP3 in western blot and immunohistochemistry.

Applications

Immunohistochemistry - This antibody can be used at 5-10 µg/ml with the appropriate secondary reagents to detect human sFRP3 on the formalin fixed, paraffin-embedded human kidney tissues (ABC).

Immunofluorescence- This antibody can be used at 5-10 µg/ml with the appropriate secondary reagents to detect human sFRP3 on methanol-fixed HeLa cells.

Western Blot-This antibody can be used at 5-10 µg/ml with the appropriate secondary reagents to detect human sFRP3 on the 293T cell lysate and HeLa whole cell lysate under reduce condition.

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Info@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.